Usage prévu

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Intended purpose

Ninonasal® is an immunochromatographic assay for the qualitative detection of SARS-CoV-2 antigen in nasal and nasopharyngeal specimens. Ninonasal® is for professional use only and is intended for use as an aid in the diagnosis of SARS-CoV-2 infection. NG Biotech determined the analytical performance characteristics of the assay based on nasopharyngeal specimens. Note that the concentration of viral antigen may vary in the course of the disease and might fall below the detection threshold of the test, particularly when the duration of disease symptoms increases (specimens collected within 7 days of the first symptom are generally more likely to give a positive result). Infectivity cannot be excluded on the basis of negative test results alone.

Summary

Early in January 2020, a new coronavirus (named "SARS-CoV-2") was identified as being the infectious agent causing a viral pneumonia outbreak in Wuhan, China, where the first cases experienced symptoms in December 2019. The disease it causes was named "coronavirus disease 2019" (in short "COVID-19"). Coronaviruses are enveloped RNA viruses widely present among human beings, other mammals, and birds, and lead to respiratory, enteric, hepatic, and neurological diseases. Six coronavirus species are known to cause human diseases. Four viruses - 229E, OC43, NL63, and HKU1 - are common and typically cause common cold symptoms in immunocompetent people. Two other strains - severe acute respiratory syndrome coronavirus (SARS-COV) and Middle East respiratory syndrome coronavirus (MERS-COV) - are of zoonotic origin and have been associated with some potentially life-threatening diseases. Common COVID-19 infection signs include respiratory symptoms, fever, cough, shortness of breath, and breathing difficulties. In more serious cases, the infection may cause pneumonia, severe acute respiratory syndrome, renal failure, and even death. The standard recommendations to prevent infection spread include regular hand washing, cover of mouth and nose when coughing and sneezing. Avoid any close contact with any person experiencing respiratory disease symptoms such as cough and sneezing.

Principle

During the test, the sample reacts with colloidal gold nanoparticles conjugated to monoclonal antibodies directed against SARS-CoV-2 nucleocapsid protein. Then, the mix migrates through the membrane by capillary action, and reacts with another monoclonal antibody directed against SARS-CoV-2 nucleocapsid protein printed at the test line (T) level. If the SARS-CoV-2 concentration present in the sample is higher than the detection limit, a coloured line appears on line T. A colour line should always be present on the control line marked "C", indicating sufficient amount was applied and the test migration occurred properly.

Reagents and materials supplied

Each kit contains:

- 5 cassettes individually packaged in aluminium pouches with desiccant
- 5 extraction tubes and 5 dropper caps*
- 1 buffer solution vial with dropper
- 5 sterile swabs* for nasal sample collection
- 1 package insert

*Due to potential shortages of ancillary medical products related to COVID-19, the nasal swabs and extraction tubes/dropper caps provided may be subject to change.

Materials required but not supplied

- Swab for nasopharyngeal sample collection
- Timer
- Personal protective equipment

Precautions

- For professional *in-vitro* diagnostic use only.
- It is recommended to use the test under a type-II microbiological safety cabinet to protect the person handling SARS-COV-2 virus.
- Do not use the test if the pouch is torn or damaged.
- If the pouch was stored between 4 and 8°C, wait at least 10 minutes so that the test reaches room temperature (15-25°C).
- The test cassette should remain in the sealed pouch until use.
- Once the aluminium pouch is open, carry out the test rapidly.
- Do not reuse the device, nor the swab nor the extraction tube.
- Handle all the samples as if they contained infectious agents. Follow the established precautions against microbiological risks throughout all the procedures, and follow the standard procedures to collect and dispose of the samples appropriately.
- Wear protective clothing such as mask, lab coats, disposable gloves, and protective eyewear when analysing the samples.

- Make sure an appropriate amount of sample is used to carry out the tests. If the sample volume is too high or too low, the results may deviate.
- The used test should be disposed of in accordance with local regulations.
- Do not eat, drink, or smoke within the area where samples and kits are handled.

Storage and stability

Store the test in the sealed aluminium pouch between 4 and 30°C. The test is stable until the expiry date printed on the pouch. The test device should remain in the sealed pouch until use. Do not freeze. Do not use after the expiry date.


Nasal specimen collection

1. Wear protective gloves, protective eyewear, and a protective mask.
2. Bring tests to room temperature prior testing.
3. Place clean extraction tube labeled with the patient ID into the tube holder (side of the box).
4. Holding the extraction buffer bottle vertically and without touching the edge of the tube, add 15 drops into the extraction tube (~450µL). In case of foam (end of the first vial of buffer) switch to second full vial.
5. Collect the nasal specimen using the described method with the swab provided into the kit.

For nasal sample collection use the swab provided in the kit, **do not use any other swabs.**

Collection method for nasal specimen :

- Insert the swab into the right nostril (approx. 2cm) and swab the inside of the nostril with rotative movements 5 times or more during approximately 15 seconds.
- Repeat the same procedure in the left nostril.

Nasopharyngeal specimen collection

1. Wear protective gloves, protective eyewear, and a protective mask.
2. Bring tests to room temperature prior testing.
3. Place clean extraction tube labeled with the patient ID into the tube holder (side of the box).
4. Holding the extraction buffer bottle vertically and without touching the edge of the tube, add 10 drops into the extraction tube (~300µL). In case of foam (end of the first vial of buffer) switch to second full vial.
5. Collect the nasopharyngeal specimen using the standard laboratory method with a validated swab.

Nasopharyngeal specimen :

- Keep the patient's head slightly tilted backwards.
- Insert the swab into the nostril and gently push it as far as possible along the floor of the nasal cavity, parallel to the palate, as described in the picture:

- Swab the epithelium of the nasal fossae with rotating movements for 15 seconds, scraping the walls of the nasopharynx to collect cells.
- Slowly remove the swab and place it in the extraction tube.

Remarks :

Swabbing can be performed in both nostrils of the same patient using the same swab, however this is not necessary if the swab is properly saturated with fluid from the first nostril sampled. Do not use samples that are contaminated with blood.

For nasopharyngeal sample collection use one of those validated swab:

- Copan ref. 503CS01;
- Noble Bio ref. NFS-1;
- Labomoderne ref. LMR1925.


Ninonasal®

Rapid test for the qualitative detection of SARS-CoV-2 antigen in human nasal and nasopharyngeal specimens. For in vitro diagnostic professional use only.

Test procedure

The following procedure is intended for fresh sample collected directly from the patient after above procedure.

- 1 / Unload the swab immediately after sampling by rotating it in the buffer 5-10 times and squeezing it to extract as much volume of sample as possible.
- 2 / Before removing the swab, squeeze it against the tube walls to release as much liquid as possible.
- 3 / Put the swab back into its original package and dispose of it according to local regulations.
- 4 / Remove the cassette from the aluminium pouch. Place the test on a clean and flat surface.
- 5 / Immediately after extraction, attach a dropper cap to the extraction tube, add 3 drops (120µL) of the mixture in the "S/R" well of the cassette, and then start the timer.
- 6 / Read the result after exactly 15 minutes. Do not read after more than 20 minutes to avoid false positive results.
- 7 / Dispose of the test and the tube in accordance with local regulations.

Results Interpretation

Negative result

Only one red line appears within the control area (C). The sample does not contain virus or the virus concentration cannot be detected by the test. This should be interpreted as a negative result.


Positive result

Two red lines appear, one in the control area (C) and one in the test area (T). The sample contains SARS-CoV-2 virus and should be interpreted as a positive result.

Note: The intensity of the red test line (T) may vary depending on SARS-CoV-2 virus level in the sample. A low-intensity line should be considered as a positive result.

Invalid result

If the control line (C) does not appear, the test result is invalid. Most often, the control line does not appear because of insufficient sample volume or an incorrect procedure. The kit may have been deteriorated. Repeat the procedure using a new test. If the problem persists, do not reuse the kit and contact your distributor.

Quality control

An internal control is included in the test. When the control line appears, it confirms the sample volume was sufficient and the test migration occurred correctly. No control sample is provided with this kit.

Limitations

- 1 / Ninonasal® is for professional in vitro diagnostic use only.
- 2 / This test does not allow to establish the quantitative SARS-CoV-2 value in the sample.
- 3 / This test allows to confirm the presence of SARS-CoV-2 nucleocapsid antigen in a symptomatic patient. In case of negative result **it shouldn't be used as the only diagnostic criterion for SARS-COV-2 infection.**
- 4 / As for any diagnostic test, the doctor should interpret the results taking into account the other clinical data.
- 5 / If the test result is negative, a confirmation test by PCR should be performed.
- 6 / The test will produce a negative result if the concentration of SARS-CoV-2 virus in the sample is lower than the limit of detection or if the nasopharyngeal or nasal specimen collection was not performed adequately.

Performances

1. Clinical evaluation of Ninonasal®

Ninonasal® has been evaluated on 518 nasal samples confirmed by RT-PCR, including 390 negative samples and 128 positive samples. Overall specificity was evaluated at **99%** with IC95% [98-100%].

The results are presented in the following table :

		RT-PCR			
		Positive		Negative	
		Ct* Value	Ct≤30	Ct≤34	-
Ninonasal®	Positive		113	125	2
	Negative		0	3	388
	Total		113	128	390
		SENSITIVITY		SPECIFICITY	
		100%	98%	99%	
		IC95% [97-100]	IC95% [93-100]	IC95% [98-100]	

*Ct = Cycle threshold

Clinical performances have also been evaluated on 117 nasopharyngeal samples confirmed by RT-PCR, including 93 negative samples and 84 positive samples collected with validated swab.

Overall specificity was evaluated at **99%** with IC95% [94-100%]. Sensitivity was evaluated at 89% with IC95% [80-95%] for samples below Ct 30 and at **96%** with IC95% [88-100%] for samples below Ct 25.

2. Limit of detection

The limit of detection was evaluated at 1ng/mL with recombinant nucleocapsid protein. Higher concentration will lead to positive results when processed on the assay.

3. Repeatability and reproducibility

Repeatability and reproducibility of Ninonasal® were evaluated with positive and negative samples. No significant difference was observed with between lots, between operators and inter series results.

4. Interferences

The following transport media have no impact on Ninonasal®: UTM (Copan, ref 346C), physiological serum, demineralized water.

5. Cross reactions

a. Human corona viruses

No cross reaction was observed with the viral nucleoprotein of the coronavirus HCoV-N63, MERS-CoV, HCoV-OC43, HCoV-229E strains up to a concentration of 1µg/mL (i.e. equivalent to about 12,04 x 10⁹ particules/mL). The viral nucleoprotein of SARS-CoV coronavirus strain cross-reacts with the Ninonasal® for concentrations superior or equal to 5ng/mL (i.e. superior or equal to about 6,02 x 10⁶ particules/mL).

b. UK Strain VOC202012/01

The results from the study conducted on the recombinant nucleocapsid from VOC202012/1 (also known as B.1.1.7 or « UK strain ») show that the monoclonal antibodies used in Ninonasal® are able to detect this variant of concern.

c. Miscellaneous

Cross reactions on Ninonasal® were evaluated by testing the following bacteria and viruses: enterovirus, respiratory syncytial virus, adenovirus, influenza B, parainfluenza, mycoplasma pneumoniae, Bordetella pertussis, human metapneumovirus, influenza A HA and influenza A H3, influenza A H1 2009. Those viruses have no impact in the results of the Ninonasal®.

Bibliography

1. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. Nat Rev Microbiol 2019;17:181-192. PMID:30531947. DOI:10.1038/s41579-018-0118-9
2. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. Trends Microbiol 2016;24:490-502. PMID:27012512. DOI:10.1016/j.tim.2016.03.003.
3. Weiss SR, Leibowitz JL. Coronavirus pathogenesis. Adv Virus Res 2011;81:85-164. PMID:22094080. DOI:10.1016/B978-0-12-385885-6.00009-2.
4. World Health Organization (WHO). WHO Statement Regarding Cluster of Pneumonia Cases in Wuhan, China. Beijing: WHO; 9 Jan 2020. [Accessed 26 Jan 2020]. <https://www.who.int/china/news/detail/09-01-2020-who-statement-regarding-cluster-of-pneumonia-cases-in-wuhan-china>.
5. World Health Organization (WHO). Coronavirus. <https://www.who.int/health-topics/coronavirus>.

Symbols

	Content sufficient for		Expiry date
	In vitro Diagnostic Medical Device		Do not reuse
	Batch number		Product reference
	Read the IFU before use		Temperatures limits
	Manufacturer		Do not use if damaged packaging
	Keep dry		

NG Biotech
Z.A. Courbouton, Secteur 1
35480 Guipry France
Phone: +33 (0)2 23 30 17 83
Fax: +33 (0)9 71 70 53 10
Email: contact@ngbiotech.com